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The present document contains data on the technical specifications of the product with characteristics, necessary for correct and safe product maintenance. Study the present manual carefully, prior to the beginning of its use or storage and/or setting-up read this manual thoroughly as directed. Provide availability of the present manual to the personnel using the given product. Examine all instructions and maintenance operating guide on equipment and materials which are used together with the present product.

## 1 SYMBOLIC ICONS AND MARKS

|   |   |
|---|---|
|    | Manufacturer  |
|    | Authorized Representative in the European Community |
|    | Handle With Care                                    |
|    | Fragile   |
|    | Keep Dry  |
|    | This Way Up   |
|    | Caution, Consult Accompanying Documents             |
|    | Date Of Manufacture                                 |
|    | Use By  |
|    | Batch Code  |
|   | Serial Number                                       |
|  | Do Not Reuse  |
|  | Sterilization using irradiation                     |
|  | Temperature Limitation                              |

## 2 PURPOSE

Single use set for active plasmapheresis (is used in combination with the device for membranous plasmapheresis AMPlD-«TT» «HEMOFENIX», manufactured by «TRACKPORE TECHNOLOGY Corporation») and set for passive plasmapheresis (further sets for passive and active plasmapheresis) are used for haemocorrection with methods of membrane plasmapheresis at acute and chronic diseases.

## 3 CONTRA-INDICATIONS FOR PLASMAPHERESIS PERFORMANCE

- Absolute contra-indications:
  - irreversible damages of a brain and other vital organs,
  - terminal conditions,
  - proceeding bleeding,
  - adequate anticoagulation can not be reached.
- Relative contra-indications:
  - high risk of internal hemorrhage,
  - heart insufficiency,

- hypoproteinemia.

## 4 RESTRAINED PLASMAPHERESIS PERFORMANCE

- Active plasmapheresis performance method:
  - Children till 1 year with weight of body up to 10 kg
- Passive plasmapheresis performance method:
  - One needle syringe plasmapheresis for newborns with body weight up to 700 gram.

## 5 BASIC SAFETY REQUIREMENTS

### WARNING

- The present product is applied only for plasmapheresis procedure performance.
- The present product is applied only on prescriptions of physicians.
- The present product is disposable.
- During medications usage, before/after product applying, the method, time and dose of their induction shall be calculated on the basis of the manufacture instruction due to the condition of the patient.
- If any deviations from normal condition with the patient during procedure are present, it is necessary to apply corresponding measures for specifying of reasons for such symptoms.
- To prevent bacterial contamination aseptic rules are strictly recommended during the product applying.

## 6 PRECAUTIONARY MEASURES

- Do not use the product after the expiry of serviceable life specified on the label.
- The careful treatment is applicable to the present product.
- Any possible impacts or product drops are not admissible.
- Strong vibrations and impacts during transportation and storage can cause the product damage.
- Do not expose the product to inner pressure higher than 100 kPa (300 mmHg).

## 7 PREPARATORY STEPS FOR PLASMAPHERESIS PROCEDURE PERFORMANCE

Applying of sets with the expired serviceable life is forbidden!

Sets for active and passive plasmapheresis shall not be used if entire consumption package is damaged!

Do not use the plasma filter with other bloodline types

Sets for active and passive plasmapheresis are due to be used strictly after the consumption package is opened.

Only specially trained physicians acknowledged with the present operation manual can have the admission to plasmapheresis procedure performance

### 7-1 List of required equipment and expendable materials for plasmapheresis procedure performance:

- AMPlD-«TT» «HEMOFENIX» device for membranous plasmapheresis;
- Set for active plasmapheresis;
- Polymeric container with isotonic of sodium chloride solution;
- Container with anticoagulant solution (ACD-A);
- Sterilizing case (sterile napkin, Bilroth's hemostatic forceps, vial with 70 % of alcohol, disposable gauze sponges and cotton balls, bandage, disposable gloves, adhesive plaster strips);
- Tonometer.


Detailed instructions on procedure are specified in the device operation manual.

### 7-2 The list of necessary equipment and expendable materials for passive plasmapheresis procedure performance:

- Set for passive plasmapheresis;
- Transfusion rack;
- All-purpose clamp;
- Polymeric container with isotonic sodium chloride solution;
- Container with anticoagulant solution (ACD-A);
- Sterilizing case (sterile napkin, Bilroth's hemostatic forceps, vial with 70% of alcohol, disposable gauze sponges and cotton balls, bandage, disposable gloves, adhesive plaster strips);
- Tonometer.

Withdrawal, return of blood and filtration are carried out only under force of gravity.

### 7-3 Preparations for passive plasmapheresis:

**Fig.1**  
**Extracorporeal passive plasmapheresis**

1-8 – Clamps, 9-11 – Filter-droppers, 12 – Infusion unit, 13 – Needle for isotonic saline intake, 14 – Needle for anticoagulant intake, 15 – Tube for blood intake, 16-19 – T-Joints, 20 – Container with isotonic saline, 21 – Container with anticoagulant, 22 – «PFM-01-TT» plasma filter, 23 – Blood container, 24 – Plasma container, 25 and 26 – Air-ducts lines, 27 – Tubing needle for plasma intake.

- (1) Prepare the working place for assembling of passive plasmapheresis procedure performance, open the case with disposable sterile items.
- (2) Check up the expiry date of plasma filter with the complete set of bloodlines.
- (3) Check the integrity of plasma filter with every single complete set of bloodlines.
- (4) Open packages, take the plasma filter with the complete set of bloodlines out and lay them out on the assembly table.
- (5) Extracorporeal circuit assembling is carried out in aseptic conditions according to the diagram on Fig.1
- (6) Fix the plasma filter on transfusion rack with the help of all-purpose clamp on the patient's level. Close all bloodlines clamps.
- (7) Remove the protective cap from the needle 13 and connect the latter with the isotonic sodium chloride solution bag 20. Put the bag into the transfusion rack's corresponding socket.
- (8) Remove the protective cap from the needle 14 and connect it with the anticoagulant bag 21. Put the latter into the transfusion rack's other socket.
- (9) Fill the filtering dropper 10 to the half, unclamp at 5 and, loosening at 3, fill the line up to the T-joint 16 with the anticoagulant. Clamp at 3.
- (10) Fill the filtering dropper 9 with isotonic sodium chloride solution to the half, loosen the clamp 2 and shift 80 ml of this solution into the blood collecting bag 23. Clamp at 2 and 5.
- (11) Remove the protective cap from the air emission branch 25 and loosen the clamp 7. Take the bag 23 by the hand with the air emission branch up and apply pressure until the air is out, then clamp at 7 and place the bag 23 at the upper part of the transfusion rack.
- (12) Fix the dropper 11 at special place on the hanger. Remove the protective cap from the air emission branch 26, loosen the clamp 8. Isotonic sodium chloride solution starts filling plasma filter and filter-dropper 11. After it is filled in up to ¼ of its volume clamp at 8.
- (13) Remove the protective cap from the infusion unit 12, unclamp at 1 and 4, and fill the rest of the line up to the infusion unit with isotonic sodium chloride solution. Clamp at 1 and 4.
- (14) (14) Remove the protective cap from the connector 27, unclamp at 6. The solution left in the bag 23 is to be used to fill the plasma compartments of PFM-01-TT with plasma efferent lines. Clamp at 6.
- (15) Replace the bag 23 to the lower position. Unclamp at 5 and, while loosening the clamp 3, shift the anticoagulant solution (ACD glucosid) from the corresponding reservoir 21 into the bag 23 in a dose volume of 80 ml. Clamp at 3 and 5.

### 8 PLASMAPHERESIS PROCEDURE PERFORMANCE:

- (1) Put a tourniquet on the patient's arm. Insert a fistular needle (catheter) into a cubital vein of the patient and attach the line's infusion bloodline unit to it. Unclamp at 1 and 5 and start collecting blood into the bag 23 (dangle from time to time to let the blood evenly mix with the anticoagulant).
- (2) When the bag 23 is full of blood, (500 ml) to untie the tourniquet, clamp at 5, unclamp at 2 and flush the infusion unit and needle with isotonic sodium chloride solution.
- (3) Then clamp at 1, unclamp at 5 and fill the rest of the line up to the bag 23 with isotonic sodium chloride solution, then clamp at 2 and 5 and move the bag 23 to the top position of the transfusion rack.
- (4) Unclamp at 1 and 4, and once the plasma filter is filled with blood, unclamp at 6 on the plasma withdrawal line 28 (during the second and following cycles unclamp at 4 and at 6 simultaneously). Blood starts to move along the plasma filter's blood compartments.
- (5) During blood recirculation and plasma separation the clamp at 2 may be somewhat loosened to let in isotonic sodium chloride solution or another plasma substitute at a speed equaling the speed of plasma separation.
- (6) After the blood in the bag 23 is fully emptied, clamp at 6 and later clamp at 4. After clamp at 2 is loosened, the infusion unit and the needle must be flushed by isotonic sodium chloride. Clamp at 1 and at 2.
- (7) Move the bag 23 down, unclamp at 5 and let in a portion of the anticoagulant by loosening the clamp 3. Clamp at 3. Put on a tourniquet, unclamp at 1 and resume blood withdrawal. Alternate withdrawal and returning cycles until the end of plasmapheresis procedure performance.
- (8) The received plasma is being collected into the container for plasma collecting 24 with removal of extra air through airway.

### 9 FINAL PROCEDURE STEP

- (1) Once the planned plasma volume is obtained and the bag 23 is emptied, move the latter down, unclamp at 2 and 5, and fill the bag with the remaining volume of isotonic sodium chloride solution or another plasma substitute. Clamp at 2 and at 5.
- (2) Move the bag 23 up. Unclamp at 1 and 4 to return the blood remaining in the system and fully replace the extracted plasma volume.
- (3) Clamp at 4 and 1, remove the needle (catheter) from the vein. Applying an aseptic bandage on the site of veinpuncture makes the procedure complete.

### WARNING

- It is very important to watch the air not to trap into the system during the plasmapheresis procedure performance, it will cause high risk of air impact into the blood flow.
- The special attention must be focused on the dose of anticoagulant considering blood clotting time of a patient.
- During the procedure of the device application some part of blood protein is being removed and that may cause hypostasis or pressure drop as a result of hypoproteinemia. Therefore, use the replacing colloid solution if needed.

### 10 SET COMPLEXES

#### 10.1 Active plasmapheresis set

| Title   | Quantity, pcs. |
|---|----------------|
| Membranous plasma filter disposable sterile «PFM-01-TT» | 1              |
| Blood line for active plasmapheresis                    | 1              |
| Individual pack unit                                    | 1              |
| Consuming package                                       | 1              |
| Corrugated box  | 1              |
| Present instructions for use                            | 1              |

#### 10.2 Passive plasmapheresis set

| Title   | Quantity, pcs. |
|---|----------------|
| Membranous plasma filter disposable sterile «PFM-01-TT» | 1              |
| Blood line for passive plasmapheresis                   | 1              |
| Individual pack unit                                    | 1              |
| Consuming package                                       | 1              |
| Corrugated box  | 1              |
| Present instructions for use                            | 1              |

### 11 SPECIFICATIONS AND TECHNICAL CHARACTERISTICS

## 11.1 Specification

### Membranous plasma filter disposable sterile PFM-01-TT

|   |  |
|---|--|
| <b>Part inside</b>                                |  |
| High-density polyethylene                         | Cap and bottom base of pressure module, tube plug  |
| Polyethylene foil                                 | Blood frame and plasma frame                       |
| Stockinet fabric from synthetic monofilament yarn | Net (net frame)                                    |
| Polyethyleneterephthalate (lavsan)                | Track membrane                                     |
| <b>Package body</b>                               |  |
| Shockproof polystyrene                            | Decorative cap, decorative bottom base. Cover plug |

### Blood lines for active and passive plasmapheresis

|   |  |
|---|--|
| Granulous medical PVC plasticate                                    | Drop chamber. Pipe with diameter 3,3, plasma collecting containers, measuring device |
| Medical soft tubing plasticate                                      | Air trap frame<br>Cap of air trap, connector for patient, Pipe with diameter 4,3     |
| Super Concentrate   | Clip   |
| Sieve material made from mixed fabrics: silky and synthetic fabrics | Fluid filter   |
| Shockproof polystyrene  | Polystyrene (polymeric) combined needle, adaptor socket. Clip frame, roller          |
| Low-density polyethylene  | Filter, T-joint, drop accumulator  |
| High -density polyethylene  | Capsule, tap, cap of patient`s connector, fixing lock. Consumption package (pack)    |

## 11.2 Specifications

| Basic function  | Active and passive plasmapheresis performance |
|---|---|
| Set weight( gr.), not more  | 700   |
| Volume of initial filling of extracorporeal circuit(ml) must not exceed   | 60  |
| Volume of initial filling of blood chambers shall not exceed  | 20  |
| Pumping length, mm  | 202   |
| Inner diameter of pumping length, mm  | 12.7  |
| Maximum pressure kPa(mmHg)  | 40 (300)                                      |
| Effective membrane area, m <sup>2</sup>   | 0,15  |
| Duration time (min)of application, shall not exceed   | 60  |
| Receiving of plasma/hour, active plasmapheresis procedure performance, at 40% of original haematocrit and 10% of blood (ml) predilution not less than                                   | 1000  |
| YReceiving of plasma/hour, passive plasmapheresis procedure performance, at 40% of original haematocrit and 10% of predilution with 1 m of height under force of gravity, not less than | 500   |
| The amount of following particles in received plasma must not exceed  |   |
| Erythrocytes  | 0,1x10 <sup>12</sup>                          |
| Leucocyte   | 0,1x10 <sup>9</sup>                           |
| Thrombocytes  | 10x10 <sup>9</sup>                            |
| The level of free haemoglobin both in the reinfused blood and removed plasma must not exceed regarding the original blood, more than  | 0,025   |
| Methods of sterilization  | Gamma radiation                               |

Technical specifications of the present manual comply with the its edition period and due to be changed.

## 12 RULES OF TRANSPORTATION AND STORAGE

Sets for active and passive plasmapheresis are transported by all kinds of covered vehicles in comply with requirements and rules of freight transportation working on transport of each type.

Active and passive sets are stored in dry premises on racks, at distance not less than 1 m from heating devices where the temperature from plus 5° C to plus 40° C is maintained, meanwhile it is necessary to avoid hit of direct solar beams and an ultraviolet.

## 13 RECYCLING INFORMATION

Products contacting the blood of the patient shall be deleted in compliance with a standard procedure, authorized for work with pathological waste in medical institution.

### WARNING

During recycling of sets for active and passive plasmapheresis after contacting the patient, special care is required to avoid personnel contamination.

## 14 MANUFACTURER WARRANTY


- «TRACKPORE TECHNOLOGY Corporation» guarantees the compliance of sets for active and passive plasmapheresis to ISO 13960 requirements and the Instruction № 93/42/EEC about medical product.
- «TRACKPORE TECHNOLOGY Corporation» undertakes to replace the sets for active and passive plasmapheresis , if any defects revealed.
- The consumer undertakes to ascertain in conformity of a product to offered applying and to evaluate all risks prior to start the appliance of sets.
- In case of the reuse applying of sets for active and passive plasmapheresis «TRACKPORE TECHNOLOGY Corporation» shall not guarantee safety and shall not bear responsibility for consequences.
- Sets for active and passive plasmapheresis are inner sterile, pyrogen-free and non-toxic.

The above-stated product of medical equipment corresponds to ISO 13960 requirements and 93/42/EEC Instructions about medical products. Has been inspected according to the Appendix II, proved by registering authority.

## 15 MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

«TRACKPORE TECHNOLOGY Corporation»  
Priborostroiteley str. 3G, building 1, Dubna, Moscow Region, Russian Federation  
Tel./Fax: +74957777957  
www.trackpore.ru  
e-mail: trackpore@trackpore.ru

PETER DEVELOPMENT COMPANY Ltd.  
Bulgaria, Burgas; 8000, Burgas District; Burgas Municipality; Lazur Municipal District,  
11 Aboba str., floor 4  
Telephone: + 35955422534  
e-mail: peter\_dk@abv.bg